

RT² Profiler PCR Array (96-Well Format and 384-Well [4 x 96] Format)

Human EGF / PDGF Signaling Pathway

Cat. no. 330231 PAHS-040ZA

For pathway expression analysis

Format	For use with the following real-time cyclers
RT ² Profiler PCR Array, Format A	Applied Biosystems [®] models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad [®] models iCycler [®] , iQ™ 5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf [®] Mastercycler [®] ep realplex models 2, 2s, 4, 4s; Stratagene [®] models Mx3005P [®] , Mx3000P [®] ; Takara TP-800
RT ² Profiler PCR Array, Format C	Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block)
RT ² Profiler PCR Array, Format D	Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon [®] , DNA Engine Opticon 2; Stratagene Mx4000 [®]
RT ² Profiler PCR Array, Format E	Applied Biosystems models 7900HT (384-well block), ViiA 7 (384-well block); Bio-Rad CFX384™
RT ² Profiler PCR Array, Format F	Roche [®] LightCycler [®] 480 (96-well block)
RT ² Profiler PCR Array, Format G	Roche LightCycler 480 (384-well block)
RT ² Profiler PCR Array, Format H	Fluidigm [®] BioMark™



Description

The Human EGF / PDGF Signaling RT² Profiler PCR Array profiles the expression of 84 genes related to EGF / PDGF-mediated signal transduction. This array includes genes that encode members of the Epidermal Growth Factor (EGF) and Platelet-Derived Growth Factor (PDGF) Signaling Pathways. The array represents growth factors and their receptors, proteins involved in signaling downstream of the growth factors, and the genes targeted by growth factor-mediated stimulation. Members of the MAP Kinase pathway, the AKT / PI3K pathway, and other downstream signaling cascades are included as well. Using real-time PCR, you can easily and reliably analyze expression of a focused panel of genes related to EGF / PDGF signaling with this array.

For further details, consult the *RT² Profiler PCR Array Handbook*.

Shipping and storage

RT² Profiler PCR Arrays in formats A, C, D, E, F, and G are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products. RT² Profiler PCR Arrays in format H are shipped on dry ice or blue ice packs.

For long term storage, keep plates at -20°C .

Note: Ensure that you have the correct RT² Profiler PCR Array format for your real-time cycler (see table above).

Note: Open the package and store the products appropriately immediately on receipt.

Array layout (96-well)

For 384-well 4 x 96 PCR arrays, genes are present in a staggered format. Refer to the *RT² Profiler PCR Array Handbook* for layout.

	1	2	3	4	5	6	7	8	9	10	11	12
A	ACTR2	AKT1	AKT2	AKT3	ARAF	ATF1	ATF2	BAD	BCAR1	BCL2	BRAF	CASP3
B	CASP9	CBL	CCND1	CHUK	COL1A1	CREB1	CSNK2A1	CSNK2B	DUSP1	DUSP6	EGF	EGFR
C	EGR1	EIF4E	ELK1	EPS8	FASLG	FN1	FOS	FOXO3	GAB1	GRB2	GSK3A	GSK3B
D	HBEGF	HRAS	IKKB	IL2	JAK1	JUN	KRAS	LTA	MAP2K1	MAP2K4	MAP2K7	MAP3K2
E	MAPK1	MAPK10	MAPK3	MAPK8	MAPK9	MKNK1	MMP7	NCK2	NFATC3	NFKB1	NRAS	NUF62
F	PDGFA	PDGFB	PDGFRA	PDPK1	PIK3CA	PIK3R1	PIK3R2	PLAT	PLCG1	PPP2CA	PRKCA	PTEN
G	RAF1	RAP1A	RASA1	RHOA	RPS6KA5	RPS6KB1	SHC1	SRC	STAT1	STAT3	STAT5A	TP53
H	ACTB	B2M	GAPDH	HPRT1	RPLP0	HGDC	RTC	RTC	RTC	PPC	PPC	PPC

Gene table: RT² Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Hs.728857	NM_005722	ACTR2	ARP2 actin-related protein 2 homolog (yeast)
A02	Hs.525622	NM_005163	AKT1	V-akt murine thymoma viral oncogene homolog 1
A03	Hs.631535	NM_001626	AKT2	V-akt murine thymoma viral oncogene homolog 2
A04	Hs.498292	NM_005465	AKT3	V-akt murine thymoma viral oncogene homolog 3 (protein kinase B, gamma)
A05	Hs.446641	NM_001654	ARAF	V-raf murine sarcoma 3611 viral oncogene homolog
A06	Hs.648565	NM_005171	ATF1	Activating transcription factor 1
A07	Hs.592510	NM_001880	ATF2	Activating transcription factor 2
A08	Hs.370254	NM_004322	BAD	BCL2-associated agonist of cell death
A09	Hs.479747	NM_014567	BCAR1	Breast cancer anti-estrogen resistance 1
A10	Hs.150749	NM_000633	BCL2	B-cell CLL/lymphoma 2
A11	Hs.550061	NM_004333	BRAF	V-raf murine sarcoma viral oncogene homolog B1
A12	Hs.141125	NM_004346	CASP3	Caspase 3, apoptosis-related cysteine peptidase
B01	Hs.329502	NM_001229	CASP9	Caspase 9, apoptosis-related cysteine peptidase
B02	Hs.504096	NM_005188	CBL	Cas-Br-M (murine) ecotropic retroviral transforming sequence
B03	Hs.523852	NM_053056	CCND1	Cyclin D1
B04	Hs.198998	NM_001278	CHUK	Conserved helix-loop-helix ubiquitous kinase
B05	Hs.172928	NM_000088	COL1A1	Collagen, type I, alpha 1
B06	Hs.516646	NM_004379	CREB1	CAMP responsive element binding protein 1
B07	Hs.644056	NM_001895	CSNK2A1	Casein kinase 2, alpha 1 polypeptide
B08	Hs.73527	NM_001320	CSNK2B	Casein kinase 2, beta polypeptide
B09	Hs.171695	NM_004417	DUSP1	Dual specificity phosphatase 1
B10	Hs.298654	NM_001946	DUSP6	Dual specificity phosphatase 6
B11	Hs.419815	NM_001963	EGF	Epidermal growth factor
B12	Hs.488293	NM_005228	EGFR	Epidermal growth factor receptor
C01	Hs.326035	NM_001964	EGR1	Early growth response 1
C02	Hs.249718	NM_001968	EIF4E	Eukaryotic translation initiation factor 4E
C03	Hs.181128	NM_005229	ELK1	ELK1, member of ETS oncogene family
C04	Hs.591160	NM_004447	EPS8	Epidermal growth factor receptor pathway substrate 8
C05	Hs.2007	NM_000639	FASLG	Fas ligand (TNF superfamily, member 6)
C06	Hs.203717	NM_002026	FN1	Fibronectin 1
C07	Hs.728789	NM_005252	FOS	FBJ murine osteosarcoma viral oncogene homolog
C08	Hs.220950	NM_001455	FOXO3	Forkhead box O3
C09	Hs.80720	NM_002039	GAB1	GRB2-associated binding protein 1
C10	Hs.444356	NM_002086	GRB2	Growth factor receptor-bound protein 2
C11	Hs.466828	NM_019884	GSK3A	Glycogen synthase kinase 3 alpha
C12	Hs.445733	NM_002093	GSK3B	Glycogen synthase kinase 3 beta
D01	Hs.799	NM_001945	HBEGF	Heparin-binding EGF-like growth factor
D02	Hs.37003	NM_005343	HRAS	V-Ha-ras Harvey rat sarcoma viral oncogene homolog
D03	Hs.597664	NM_001556	IKKB	Inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase beta
D04	Hs.89679	NM_000586	IL2	Interleukin 2
D05	Hs.207538	NM_002227	JAK1	Janus kinase 1
D06	Hs.714791	NM_002228	JUN	Jun proto-oncogene
D07	Hs.505033	NM_004985	KRAS	V-Ki-ras2 Kirsten rat sarcoma viral oncogene homolog
D08	Hs.36	NM_000595	LTA	Lymphotoxin alpha (TNF superfamily, member 1)
D09	Hs.145442	NM_002755	MAP2K1	Mitogen-activated protein kinase kinase 1

Position	UniGene	GenBank	Symbol	Description
D10	Hs.514681	NM_003010	MAP2K4	Mitogen-activated protein kinase kinase 4
D11	Hs.531754	NM_145185	MAP2K7	Mitogen-activated protein kinase kinase 7
D12	Hs.145605	NM_006609	MAP3K2	Mitogen-activated protein kinase kinase kinase 2
E01	Hs.431850	NM_002745	MAPK1	Mitogen-activated protein kinase 1
E02	Hs.125503	NM_002753	MAPK10	Mitogen-activated protein kinase 10
E03	Hs.861	NM_002746	MAPK3	Mitogen-activated protein kinase 3
E04	Hs.138211	NM_002750	MAPK8	Mitogen-activated protein kinase 8
E05	Hs.484371	NM_002752	MAPK9	Mitogen-activated protein kinase 9
E06	Hs.371594	NM_003684	MKNK1	MAP kinase interacting serine/threonine kinase 1
E07	Hs.2256	NM_002423	MMP7	Matrix metalloproteinase 7 (matrilysin, uterine)
E08	Hs.529244	NM_003581	NCK2	NCK adaptor protein 2
E09	Hs.632209	NM_004555	NFATC3	Nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 3
E10	Hs.654408	NM_003998	NFKB1	Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
E11	Hs.486502	NM_002524	NRAS	Neuroblastoma RAS viral (v-ras) oncogene homolog
E12	Hs.467133	NM_012346	NUP62	Nucleoporin 62kDa
F01	Hs.535898	NM_002607	PDGFA	Platelet-derived growth factor alpha polypeptide
F02	Hs.1976	NM_002608	PDGFB	Platelet-derived growth factor beta polypeptide
F03	Hs.74615	NM_006206	PDGFRA	Platelet-derived growth factor receptor, alpha polypeptide
F04	Hs.459691	NM_002613	PDPK1	3-phosphoinositide dependent protein kinase-1
F05	Hs.553498	NM_006218	PIK3CA	Phosphoinositide-3-kinase, catalytic, alpha polypeptide
F06	Hs.132225	NM_181504	PIK3R1	Phosphoinositide-3-kinase, regulatory subunit 1 (alpha)
F07	Hs.371344	NM_005027	PIK3R2	Phosphoinositide-3-kinase, regulatory subunit 2 (beta)
F08	Hs.491582	NM_000930	PLAT	Plasminogen activator, tissue
F09	Hs.268177	NM_002660	PLCG1	Phospholipase C, gamma 1
F10	Hs.483408	NM_002715	PPP2CA	Protein phosphatase 2, catalytic subunit, alpha isozyme
F11	Hs.531704	NM_002737	PRKA	Protein kinase C, alpha
F12	Hs.500466	NM_000314	PTEN	Phosphatase and tensin homolog
G01	Hs.159130	NM_002880	RAF1	V-raf-1 murine leukemia viral oncogene homolog 1
G02	Hs.190334	NM_002884	RAP1A	RAP1A, member of RAS oncogene family
G03	Hs.664080	NM_002890	RASA1	RAS p21 protein activator (GTPase activating protein) 1
G04	Hs.247077	NM_001664	RHOA	Ras homolog gene family, member A
G05	Hs.510225	NM_004755	RPS6KA5	Ribosomal protein S6 kinase, 90kDa, polypeptide 5
G06	Hs.463642	NM_003161	RPS6KB1	Ribosomal protein S6 kinase, 70kDa, polypeptide 1
G07	Hs.433795	NM_003029	SHC1	SHC (Src homology 2 domain containing) transforming protein 1
G08	Hs.195659	NM_005417	SRC	V-src sarcoma (Schmidt-Ruppin A-2) viral oncogene homolog (avian)
G09	Hs.642990	NM_007315	STAT1	Signal transducer and activator of transcription 1, 91kDa
G10	Hs.463059	NM_003150	STAT3	Signal transducer and activator of transcription 3 (acute-phase response factor)
G11	Hs.437058	NM_003152	STAT5A	Signal transducer and activator of transcription 5A
G12	Hs.654481	NM_000546	TP53	Tumor protein p53
H01	Hs.520640	NM_001101	ACTB	Actin, beta
H02	Hs.534255	NM_004048	B2M	Beta-2-microglobulin
H03	Hs.592355	NM_002046	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
H04	Hs.412707	NM_000194	HPRT1	Hypoxanthine phosphoribosyltransferase 1
H05	Hs.546285	NM_001002	RPLP0	Ribosomal protein, large, P0
H06	N/A	SA_00105	HGDC	Human Genomic DNA Contamination
H07	N/A	SA_00104	RTC	Reverse Transcription Control
H08	N/A	SA_00104	RTC	Reverse Transcription Control
H09	N/A	SA_00104	RTC	Reverse Transcription Control
H10	N/A	SA_00103	PPC	Positive PCR Control
H11	N/A	SA_00103	PPC	Positive PCR Control
H12	N/A	SA_00103	PPC	Positive PCR Control

Related products

For optimal performance, RT² Profiler PCR Arrays should be used together with the RT² First Strand Kit for cDNA synthesis and RT² SYBR[®] Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT ² First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT ² SYBR Green qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with real-time cyclers that do not require a reference dye, including: Bio-Rad models CFX96, CFX384, DNA Engine Opticon 2; Bio-Rad/MJ Research Chromo4; Roche LightCycler 480 (96-well and 384-well); all other cyclers	330500
RT ² SYBR Green ROX™ qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Applied Biosystems models 5700, 7000, 7300, 7500 [Standard and FAST], 7700, 7900HT 96-well block [Standard and FAST] and 384-well block, StepOnePlus; Eppendorf Mastercycler ep realplex models 2, 2S, 4, 4S; Stratagene models Mx3000P, Mx3005P, Mx4000; Takara TP-800	330520
RT ² SYBR Green Fluor qPCR Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the following real-time cyclers: Bio-Rad models iCycler, iQ5, MyiQ, MyiQ2	330510

* Larger kit sizes available; please inquire.

RT² Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN® (QIAGEN Group); Applied Biosystems®, ViiA™, StepOnePlus™, ROX™ (Applied Biosystems Corporation or its subsidiaries); Bio-Rad®, iCycler®, iQ™, MyiQ™, Chromo4™, CFX96™, DNA Engine Opticon®, CFX384™ (Bio-Rad Laboratories, Inc.); Stratagene®, Mx3005P®, Mx3000P®, Mx4000® (Stratagene); Eppendorf®, Mastercycler® (Eppendorf AG); Roche®, LightCycler® (Roche Group); Fluidigm® BioMark™ (Fluidigm Corporation); SYBR® (Molecular Probes, Inc.).

1066029 03/2011 © 2011 QIAGEN, all rights reserved.

www.qiagen.com

Canada ■ 800-572-9613

Ireland ■ 1800 555 049

Norway ■ 800-18859

China ■ 8621-3865-3865

Italy ■ 800-787980

Singapore ■ 1800-742-4368

Denmark ■ 80-885945

Japan ■ 03-6890-7300

Spain ■ 91-630-7050

Australia ■ 1-800-243-800

Finland ■ 0800-914416

Korea (South) ■ 080-000-7145

Sweden ■ 020-790282

Austria ■ 0800/281010

France ■ 01-60-920-930

Luxembourg ■ 8002 2076

Switzerland ■ 055-254-22-11

Belgium ■ 0800-79612

Germany ■ 02103-29-12000

Mexico ■ 01-800-7742-436

UK ■ 01293-422-911

Brazil ■ 0800-557779

Hong Kong ■ 800 933 965

The Netherlands ■ 0800 0229592

USA ■ 800-426-8157



Sample & Assay Technologies