

QuantiNova® LNA® Probe PCR Focus Panels (96-Well Format and 384-Well [4 x 96] Format)

Human Atherosclerosis

Cat. no. 249955 UPHS-038ZA

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA Probe PCR Focus Panels are shipped at room temperature. Immediately upon receipt, they should be stored protected from light at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA Probe PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova Probe PCR Kit (Mastermix) for PCR.

Panel layout (96-well): QuantiNova LNA Probe PCR Focus Panel

For the 384-well (4 × 96) PCR panels, genes are present in a staggered format. Refer to the QuantiNova LNA Probe PCR Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	ABCA1	ACE	APOA1	APOB	APOE	BAX	BCL2	BCL2A1	BCL2L1	BID	BIRC3	CCL2
B	CCL5	CCR1	CCR2	CD44	CDH5	CFLAR	COL3A1	CSF1	CSF2	CCN2	EGR1	ELN
C	ENG	FABP3	FAS	FGA	FGF2	FN1	HBEGF	ICAM1	IFNAR2	IFNG	IL1A	IL1R1
D	IL1R2	IL2	IL3	IL4	IL5	ITGA2	ITGA5	ITGAX	ITGB2	KDR	KLF2	LAMA1
E	LDLR	LIF	LPA	LPL	MMP1	MMP3	MSR1	NFKB1	NOS3	NPY	NR1H3	PDGFA
F	PDGFB	PDGFRB	PLIN2	PPARA	PPARD	PPARG	PTGS1	RXRA	SELE	SELL	SEPLG	SERPINB2
G	SERPINE1	SOD1	SPP1	TGFB1	TGFB2	THBS4	TNC	TNF	TNFAIP3	VCAM1	VEGFA	VWF
H	ACTB	B2M	GAPDH	HPRT1	RPLP0	HGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA Probe PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	UPFH0384523	ENST00000423487.6	ABCA1	ENSG00000165029	ATP binding cassette subfamily A member 1 Source HGNC Symbol Acc HGNC 29
A02	UPFH1132210	ENST00000428043.5	ACE	ENSG00000159640	angiotensin I converting enzyme Source HGNC Symbol Acc HGNC 2707
A03	UPFH0042718	ENST00000375329.6	APOA1	ENSG00000118137	apolipoprotein A1 Source HGNC Symbol Acc HGNC 600
A04	UPFH0387821	ENST00000399256.4	APOB	ENSG00000084674	apolipoprotein B Source HGNC Symbol Acc HGNC 603
A05	UPFH0539669	ENST00000425718.1	APOE	ENSG00000130203	apolipoprotein E Source HGNC Symbol Acc HGNC 613
A06	UPFH0540159	ENST00000293288.12	BAX	ENSG00000087088	BCL2 associated X, apoptosis regulator Source HGNC Symbol Acc HGNC 959
A07	UPFH1132900	ENST00000333681.5	BCL2	ENSG00000171791	BCL2, apoptosis regulator Source HGNC Symbol Acc HGNC 990
A08	UPFH1132270	ENST00000335661.6	BCL2A1	ENSG00000140379	BCL2 related protein A1 Source HGNC Symbol Acc HGNC 991
A09	UPFH1132271	ENST00000376062.6	BCL2L1	ENSG00000171552	BCL2 like 1 Source HGNC Symbol Acc HGNC 992
A10	UPFH0184066	ENST00000622694.4	BID	ENSG00000015475	BH3 interacting domain death agonist Source HGNC Symbol Acc HGNC 1050
A11	UPFH0113281	ENST00000615299.4	BIRC3	ENSG00000023445	baculoviral IAP repeat containing 3 Source HGNC Symbol Acc HGNC 591
A12	UPFH1132783	ENST00000225831.4	CCL2	ENSG00000108691	C-C motif chemokine ligand 2 Source HGNC Symbol Acc HGNC 10618
B01	UPFH1132786	ENST00000603197.6	CCL5	ENSG000000271503	C-C motif chemokine ligand 5 Source HGNC Symbol Acc HGNC 10632
B02	UPFH0327828	ENST00000296140.4	CCR1	ENSG00000163823	C-C motif chemokine receptor 1 Source HGNC Symbol Acc HGNC 1602
B03	UPFH0175349	ENST00000445132.2	CCR2	ENSG00000121807	C-C motif chemokine receptor 2 Source HGNC Symbol Acc HGNC 1603
B04	UPFH0253499	ENST00000428726.7	CD44	ENSG00000026508	CD44 molecule (Indian blood group) Source HGNC Symbol Acc HGNC 1681
B05	UPFH1132306	ENST00000563425.2	CDH5	ENSG00000179776	cadherin 5 Source HGNC Symbol Acc HGNC 1764
B06	UPFH1132312	ENST00000479953.6	CFLAR	ENSG00000003402	CASP8 and FADD like apoptosis regulator Source HGNC Symbol Acc HGNC 1876
B07	UPFH0067087	ENST00000304636.7	COL3A1	ENSG00000168542	collagen type III alpha 1 chain Source HGNC Symbol Acc HGNC 2201
B08	UPFH1132338	ENST00000329608.11	CSF1	ENSG00000184371	colony stimulating factor 1 Source HGNC Symbol Acc HGNC 2432
B09	UPFH1132793	ENST00000296871.4	CSF2	ENSG00000164400	colony stimulating factor 2 Source HGNC Symbol Acc HGNC 2434
B10	UPFH1132340	ENST00000367976.4	CCN2	ENSG00000118523	cellular communication network factor 2 Source HGNC Symbol Acc HGNC 2500
B11	UPFH0558832	ENST00000239938.5	EGR1	ENSG00000120738	early growth response 1 Source HGNC Symbol Acc HGNC 3238
B12	UPFH1172924	ENST00000380553.8	ELN	ENSG00000049540	elastin Source HGNC Symbol Acc HGNC 3327
C01	UPFH0535657	ENST00000344849.4	ENG	ENSG00000106991	endoglin Source HGNC Symbol Acc HGNC 3349
C02	UPFH0313060	ENST00000373713.6	FABP3	ENSG00000121769	fatty acid binding protein 3 Source HGNC Symbol Acc HGNC 3557
C03	UPFH1132395	ENST00000357339.6	FAS	ENSG00000026103	Fas cell surface death receptor Source HGNC Symbol Acc HGNC 11920
C04	UPFH0337891	ENST00000622532.1	FGA	ENSG00000171560	fibrinogen alpha chain Source HGNC Symbol Acc HGNC 3661
C05	UPFH0613093	ENST00000264498.7	FGF2	ENSG00000138685	fibroblast growth factor 2 Source HGNC Symbol Acc HGNC 3676
C06	UPFH0605066	ENST00000336916.8	FN1	ENSG00000115414	fibronectin 1 Source HGNC Symbol Acc HGNC 3778
C07	UPFH0129784	ENST00000230990.7	HBEGF	ENSG00000113070	heparin binding EGF like growth factor Source HGNC Symbol Acc HGNC 3059
C08	UPFH1132462	ENST00000264832.8	ICAM1	ENSG00000090339	intercellular adhesion molecule 1 Source HGNC Symbol Acc HGNC 5344
C09	UPFH0538480	ENST00000382264.7	IFNAR2	ENSG00000159110	interferon alpha and beta receptor subunit 2 Source HGNC Symbol Acc HGNC 5433
C10	UPFH1132473	ENST00000229135.4	IFNG	ENSG00000111537	interferon gamma Source HGNC Symbol Acc HGNC 5438
		ENST00000263		ENSG000000	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	UPFH0436255	339.3	IL1A	115008	interleukin 1 alpha Source HGNC Symbol Acc HGNC 5991
C12	UPFH1132482	ENST00000442590.5	IL1R1	ENSG00000115594	interleukin 1 receptor type 1 Source HGNC Symbol Acc HGNC 5993
D01	UPFH0176161	ENST00000441002.1	IL1R2	ENSG00000115590	interleukin 1 receptor type 2 Source HGNC Symbol Acc HGNC 5994
D02	UPFH0116492	ENST00000226730.4	IL2	ENSG00000109471	interleukin 2 Source HGNC Symbol Acc HGNC 6001
D03	UPFH0282899	ENST00000296870.2	IL3	ENSG00000164399	interleukin 3 Source HGNC Symbol Acc HGNC 6011
D04	UPFH0226437	ENST00000231449.7	IL4	ENSG00000113520	interleukin 4 Source HGNC Symbol Acc HGNC 6014
D05	UPFH1132811	ENST00000231454.6	IL5	ENSG00000113525	interleukin 5 Source HGNC Symbol Acc HGNC 6016
D06	UPFH1132496	ENST00000296585.10	ITGA2	ENSG00000164171	integrin subunit alpha 2 Source HGNC Symbol Acc HGNC 6137
D07	UPFH1132497	ENST00000293379.9	ITGA5	ENSG00000161638	integrin subunit alpha 5 Source HGNC Symbol Acc HGNC 6141
D08	UPFH0136150	ENST00000562522.2	ITGAX	ENSG00000140678	integrin subunit alpha X Source HGNC Symbol Acc HGNC 6152
D09	UPFH1132499	ENST00000397850.6	ITGB2	ENSG00000160255	integrin subunit beta 2 Source HGNC Symbol Acc HGNC 6155
D10	UPFH0596732	ENST00000263923.5	KDR	ENSG00000128052	kinase insert domain receptor Source HGNC Symbol Acc HGNC 6307
D11	UPFH0347180	ENST00000592003.1	KLF2	ENSG00000127528	Kruppel like factor 2 Source HGNC Symbol Acc HGNC 6347
D12	UPFH1132819	ENST00000389658.4	LAMA1	ENSG00000101680	laminin subunit alpha 1 Source HGNC Symbol Acc HGNC 6481
E01	UPFH0616676	ENST00000252444.9	LDLR	ENSG00000130164	low density lipoprotein receptor Source HGNC Symbol Acc HGNC 6547
E02	UPFH1132822	ENST00000249075.4	LIF	ENSG00000128342	LIF, interleukin 6 family cytokine Source HGNC Symbol Acc HGNC 6596
E03	UPFH0349322	ENST00000447678.2	LPA	ENSG00000198670	lipoprotein(a) Source HGNC Symbol Acc HGNC 6667
E04	UPFH1132523	ENST00000311322.10	LPL	ENSG00000175445	lipoprotein lipase Source HGNC Symbol Acc HGNC 6677
E05	UPFH0322484	ENST00000315274.7	MMP1	ENSG00000196611	matrix metalloproteinase 1 Source HGNC Symbol Acc HGNC 7155
E06	UPFH1132827	ENST00000299855.10	MMP3	ENSG00000149968	matrix metalloproteinase 3 Source HGNC Symbol Acc HGNC 7173
E07	UPFH0191268	ENST00000262101.10	MSR1	ENSG00000038945	macrophage scavenger receptor 1 Source HGNC Symbol Acc HGNC 7376
E08	UPFH1132828	ENST00000226574.9	NFKB1	ENSG00000109320	nuclear factor kappa B subunit 1 Source HGNC Symbol Acc HGNC 7794
E09	UPFH1132897	ENST00000297494.8	NOS3	ENSG00000164867	nitric oxide synthase 3 Source HGNC Symbol Acc HGNC 7876
E10	UPFH0380153	ENST00000407573.5	NPY	ENSG00000122585	neuropeptide Y Source HGNC Symbol Acc HGNC 7955
E11	UPFH0556560	ENST00000616973.4	NR1H3	ENSG00000025434	nuclear receptor subfamily 1 group H member 3 Source HGNC Symbol Acc HGNC 7966
E12	UPFH1132608	ENST00000354513.9	PDGFA	ENSG00000197461	platelet derived growth factor subunit A Source HGNC Symbol Acc HGNC 8799
F01	UPFH0524988	ENST00000381551.8	PDGFB	ENSG00000100311	platelet derived growth factor subunit B Source HGNC Symbol Acc HGNC 8800
F02	UPFH1132610	ENST00000261799.9	PDGFRB	ENSG00000113721	platelet derived growth factor receptor beta Source HGNC Symbol Acc HGNC 8804
F03	UPFH0317135	ENST00000494753.1	PLIN2	ENSG00000147872	perilipin 2 Source HGNC Symbol Acc HGNC 248
F04	UPFH0327373	ENST00000262735.9	PPARA	ENSG00000186951	peroxisome proliferator activated receptor alpha Source HGNC Symbol Acc HGNC 9232
F05	UPFH1132629	ENST00000448077.6	PPARD	ENSG00000112033	peroxisome proliferator activated receptor delta Source HGNC Symbol Acc HGNC 9235
F06	UPFH0284890	ENST00000477039.5	PPARG	ENSG00000132170	peroxisome proliferator activated receptor gamma Source HGNC Symbol Acc HGNC 9236
F07	UPFH0450481	ENST00000619306.5	PTGS1	ENSG00000095303	prostaglandin-endoperoxide synthase 1 Source HGNC Symbol Acc HGNC 9604
F08	UPFH0039313	ENST00000356384.4	RXRA	ENSG00000186350	retinoid X receptor alpha Source HGNC Symbol Acc HGNC 10477
F09	UPFH1132838	ENST00000367777.5	SELE	ENSG00000007908	selectin E Source HGNC Symbol Acc HGNC 10718
F10	UPFH0444664	ENST00000236147.5	SELL	ENSG00000188404	selectin L Source HGNC Symbol Acc HGNC 10720

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	UPFH0032702	ENST00000550948.1	SELPLG	ENSG00000110876	selectin P ligand Source HGNC Symbol Acc HGNC 10722
F12	UPFH1132670	ENST00000299502.9	SERPINB2	ENSG00000197632	serpin family B member 2 Source HGNC Symbol Acc HGNC 8584
G01	UPFH0384736	ENST00000223095.4	SERPINE1	ENSG00000106366	serpin family E member 1 Source HGNC Symbol Acc HGNC 8583
G02	UPFH0220920	ENST00000270142.10	SOD1	ENSG00000142168	superoxide dismutase 1 Source HGNC Symbol Acc HGNC 11179
G03	UPFH0044238	ENST00000237623.11	SPP1	ENSG00000118785	secreted phosphoprotein 1 Source HGNC Symbol Acc HGNC 11255
G04	UPFH0193430	ENST00000221930.5	TGFB1	ENSG00000105329	transforming growth factor beta 1 Source NCBI gene Acc 7040
G05	UPFH1132846	ENST00000366929.4	TGFB2	ENSG00000092969	transforming growth factor beta 2 Source HGNC Symbol Acc HGNC 11768
G06	UPFH0616970	ENST00000350881.6	THBS4	ENSG00000113296	thrombospondin 4 Source HGNC Symbol Acc HGNC 11788
G07	UPFH1132731	ENST00000350763.9	TNC	ENSG00000041982	tenascin C Source HGNC Symbol Acc HGNC 5318
G08	UPFH1132978	ENST00000449264.3	TNF	ENSG00000232810	tumor necrosis factor Source HGNC Symbol Acc HGNC 11892
G09	UPFH0576817	ENST00000485192.1	TNFAIP3	ENSG00000118503	TNF alpha induced protein 3 Source HGNC Symbol Acc HGNC 11896
G10	UPFH1132856	ENST00000294728.7	VCAM1	ENSG00000162692	vascular cell adhesion molecule 1 Source HGNC Symbol Acc HGNC 12663
G11	UPFH0281656	ENST00000425836.6	VEGFA	ENSG00000112715	vascular endothelial growth factor A Source HGNC Symbol Acc HGNC 12680
G12	UPFH0533670	ENST00000261405.9	VWF	ENSG00000110799	von Willebrand factor Source HGNC Symbol Acc HGNC 12726
H01	UPFH1132936	ENST00000646664.1	ACTB	ENSG00000075624	actin beta Source HGNC Symbol Acc HGNC 132
H02	UPFH1132937	ENST00000544417.5	B2M	ENSG00000166710	beta-2-microglobulin Source HGNC Symbol Acc HGNC 914
H03	UPFH1132938	ENST00000229239.10	GAPDH	ENSG00000111640	glyceraldehyde-3-phosphate dehydrogenase Source HGNC Symbol Acc HGNC 4141
H04	UPFH1132939	ENST00000298556.8	HPRT1	ENSG00000165704	hypoxanthine phosphoribosyltransferase 1 Source HGNC Symbol Acc HGNC 5157
H05	UPFH1132941	ENST00000392514.9	RPLP0	ENSG00000089157	ribosomal protein lateral stalk subunit P0 Source HGNC Symbol Acc HGNC 10371
H06	UPFH1126608	UPL_HGDC	HGDC	UPL_HGDC	Human Genomic DNA Contamination
H07	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H08	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H09	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H10	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H11	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H12	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control



Related products

Product	Contents	Cat. no.
QuantiNova LNA Probe PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA Probe PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249945
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 μ l reactions: 20 μ l 8x gDNA Removal Mix, 10 μ l Reverse Transcription Enzyme, 40 μ l Reverse Transcription Mix (containing RT primers), 20 μ l Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova Probe RT-PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml QuantiNova Probe RT-PCR Master Mix, 20 μ l QuantiNova Probe RT Mix, 20 μ l Internal Control RNA, 500 μ l Yellow Template Dilution Buffer, 250 μ l ROX Reference Dye, 1.9 μ l RNase-Free Water	208352
QuantiNova Probe PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml 2x QuantiNova Probe PCR Master Mix, 500 μ l QuantiNova Yellow Template Dilution Buffer, 250 μ l QN ROX Reference Dye, 1.9 ml Water	208252

*Larger kit sizes available.

The QuantiNova LNA Probe PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.